

TUNEL assay using ApopTag Kit from Intergen Company (Cat #S7110)

1. **Deparaffinize tissue sections** (in a coplin jar)
 - a. Wash slides with xylene, 5 min twice.
 - b. Rehydrate in a decreasing ethanol series: 100% x2, 95% x2, 70%, 2 min each.
 - c. Rinse in PBS, 3 min x 3
2. **Pretreat tissues**
 - a. Incubate slides in Proteinase K (20 µg/ml final concentration) for 15 min at RT. To prepare 1 ml of ProK solution (good for two slides), add 4 µl 25 mg/ml ProK to 1 ml of PBS and mix. Optimal ProK treatments need to be tested for frozen sections.
 - b. Rinse well in PBS 2 x 2 min.
3. **Apply Equilibration Buffer** (blue cap)
 - a. Gently tap off excess liquid and carefully blot or aspirate around the section.
 - b. Immediately apply 75 µl Equilibration Buffer directly on the specimen, and coverslip.
 - c. Incubate for at least 10 sec at RT.
4. **Apply TdT enzyme** (for 2 slides, 33 µl TdT Enzyme + 77 µl Reaction Buffer, mix and keep on ice, good for 6hr)
 - a. Gently tap off excess liquid and carefully blot or aspirate around the section.
 - b. Immediately pipette onto the section 55 µl/5 cm² of working strength TdT enzyme and coverslip.
 - c. Incubate in a humidified chamber at 37°C for 1 hour.

5. **Apply Strop/Wash Buffer** (1ml Stop/Wash Buffer + 49 ml H₂O)
 - a. Put the slides in working strength Stop/Wash Buffer in a coplin jar), agitate for 15 sec, and incubate for 10 min at RT.
 - b. Remove an aliquot of Anti-Digoxigenin Conjugate from the stock vial sufficient to process the desired number of slides. (Warm the aliquot to RT while avoiding exposure to light.
6. **Apply Anti-Digoxigenin Conjugate** (for 2 slides, 68 μ l blocking Solution + 62 μ l Anti-Digoxigenin Conjugate, vortex and keep on ice in dark)
 - a. Wash the slides in PBS 3 x 1 min.
 - b. Gently tap off excess liquid and carefully blot or aspirate around the section.
 - c. Apply warmed (RT) working strength Anti-Digoxigenin Conjugate (rhodamin) to the slide; use about 65 μ l/5 cm² of specimen surface area and coverslip.
 - d. Incubate in a humidified chamber **in dark** for 30 min. at RT.
7. **Counter staining**
 - a. Wash slides in PBS, 2 min x 4.
 - b. Add 0.5 ml of Hoechst solution to each slides and incubate for 10 min in dark. To prepare Hoechst solution, mix 1 μ l of stock solution (10 mg/ml, Hoechst 33342, Molecular Probes, Cat. # H-3570) with TBS (PBS containing 01.% Triton-100).
 - c. Wash with PBS 5 min x 3.
 - d. Mount slides with Gel mount or ProLong Antifade Kit (Molecular Probes, Cat. # P-7481).
 - e. Seal slide with nail polisher next morning.